

## ELISA Protocol (Antibody Capture)

### Materials:

- 96-Well Microtiter Plates
- Eppendorf Tubes
- Twelve-Channel Pipettor
- 1mL Adjustable Pipettor
- Humid Chamber
- Wash Bottle or ELISA Plate Washer
- ELISA Plate Reader
- TMB

### Buffer Formulations:

#### Carbonate Buffer

*1 Liter – pH 9.5 – Store at 4°C*

- 1.59g Na<sub>2</sub>CO<sub>3</sub>
- 2.93g NaHCO<sub>3</sub>
- 2mL 10% NaN<sub>3</sub>
- DI H<sub>2</sub>O to 1L

#### Wash Buffer

- PBS
- 0.02% Thimerosal
- 0.05% Tween-20

#### Blocking Buffer

- PBS
- 0.1% BSA
- 0.02% Thimerosal
- Store in a dark bottle

### Plate Design

Mark each column with the project number. The first three rows are reserved for 3 dilutions of the pre-bleed and the following 5 rows are for dilutions of the bleed. See chart below:

	1	2	3	4	5	6	7	8	9	10	11	12		
A													1:1,000	<b>Pre-Bleed</b>
B													1:5,000	
C													1:25,000	
D													1:1,000	<b>Bleed</b>
E													1:5,000	
F													1:25,000	
G													1:125,000	
H													1:625,000	

## **DAY ONE**

### **Antigen Preparation**

1. The working peptide antigen concentration is 10µg/mL: Dissolve 1mg peptide in 200µL DI H<sub>2</sub>O (5mg/mL). Add 2µL peptide to 1mL of 50mM carbonate buffer.
2. The working protein antigen concentration is 40µg/mL: Concentrate protein to 4mg/mL and dialyze against 50mM carbonate buffer overnight. Then dilute the protein to 40ug/mL with carbonate buffer.

### **Coating**

3. Add 100µL of antigen at the recommended dilutions above to each well on the plate. Incubate at room temperature for 24 hours in a humid chamber.
  - a. Humid chamber: Place wet paper towels in the bottom of a bin and stack the plates on top of the towels. Cover with a lid.

## **DAY TWO**

### **Blocking**

4. Remove coating solution and rinse twice (2) with DI H<sub>2</sub>O
5. Flip dry on a stack of paper towels – slap to remove all traces of liquid from the wells.
6. Add 280µL blocking buffer to each well. Incubate in a humid chamber at room temperature for two (2) hours.
  - a. More blocking buffer is better, about 200-300µL per well, but make sure not to overfill each well.

### **Serum Sample**

7. Dilute serum (pre-bleed and bleed in separate Eppendorf tubes) with blocking buffer to 1:1000 (1µL serum to 1mL blocking buffer).
8. Add 125µL of 1:1000 diluted pre-immune serum to row A of the plate.
9. Using a twelve-channel pipettor, take 25µL from row A and add to row B, mixing with the pipettor. Repeat with rows B and C, removing 25µL from row B and mixing it with row C. Remove 25µL from row C and discard.
10. Add 125µL of 1:1000 diluted serum to row D of the plate.

11. Using a twelve-channel pipettor make a series of 5 dilutions for rows E through H as described above for the pre-bleed, discarding 25 $\mu$ L after row H.
12. Incubate for one hour at room temperature.

### **Secondary Antibody**

13. Dispose of the serum samples and wash the wells with the wash buffer using a wash bottle or an ELISA plate washer one (1) time. Follow with one wash using DI H<sub>2</sub>O. Drain the remaining liquid on a stack of paper towels.
14. Dilute the HRP conjugated secondary to the dilution recommended below with the blocking buffer.
  - a. Goat anti-rabbit IgG: 1:20000
  - b. anti-chicken IgY: 1:5000
  - c. anti-goat IgG: 1:10000
  - d. anti-mouse IgG: 1:5000
  - e. *Note:* Dilutions are only recommendations. Optimal dilutions will vary for each secondary and should be determined for each new lot.
15. Incubate at room temperature for one (1) hour.

### **Development**

16. Dispose of the secondary and wash the wells with the wash buffer one (1) time. Follow with one wash using DI H<sub>2</sub>O. Drain the remaining liquid on a stack of paper towels.
17. Add TMB per manufacturer instructions.
  - a. Positive wells will change to a blue color depending on signal intensity.
18. Add 90 $\mu$ L of 1M HCl to each well (or recommended media) to stop the color development and read immediately on a microtiter plate reader at A<sub>450</sub>.
  - a. The blue will change to yellow when the reaction is stopped.
19. Read multiple plates in the order the color development was stopped.
  - a. It is not recommended to have more than two plates developing at the same time. Have two timers so that the reaction is stopped after 5 minutes.