



ProSci Incorporated
12170 Flint Place
Poway, CA 92064

Toll Free: +1 (888) 513 9525
Local: +1 (858) 513 2638
Fax: +1 (858) 513 2692

techsupport@prosci-inc.com
prosci-inc.com

Breast Tissue Slide (Fibroadenoma)

CATALOG NUMBER: 10-004

Specifications

TISSUE CELL TYPE: Breast

SPECIES: Human

GENERAL DIAGNOSIS: Diseased

AVAILABLE LOTS: Lots CA1, XA1, XB1 - Fibroadenoma - Benign

Lots CB1, CC1, XC1, XD1, XF1, XG1, XH1, XJ1, XK1, XM1, XN1 - Fibroadenoma

Lot XE1 - Fibroadenoma - Tubular Adenoma

Lot XL1 - Probable Fibroadenoma - Possible Non Menstrating age

TESTED APPLICATIONS: IHC

APPLICATIONS: ProSci's human tissue slides can be used with any antibody shown to be effective in formalin-fixed tissue using a standard immunostaining protocol.

In order to avoid false negative results from DNA or RNA sequences not being present in a 4 um cellular cross-section, 10 um tissue sections are recommended. 10 um single tissue slides are available upon request. This product has not been evaluated for effectiveness in ISH or FISH applications.

Properties

FIXATION: Formalin-fixed, paraffin-embedded

STORAGE CONDITIONS: Stable for up to two years when stored in a cool, dry lab area.

Background

Please Note: ProSci now has over 700 different lots available for a wide range of tissues. Many of these lots are not listed online. If you don't see what you are looking for please call us.

Formalin-fixed, paraffin-embedded tissue is useful for a variety of applications including the targeting of proteins, RNA, and DNA by means of Immunohistochemistry (IHC), in-situ hybridization (ISH) and fluorescence in-situ hybridization (FISH), respectively. Some tissue types may require extra attention as they may contain high levels of endogenous biotin, peroxidase, or autofluorescence.

Please note: Slides shipped for your order will not be stained and are prepared for each individual order to insure quality. A 3-4 day lead time is required.

FOR RESEARCH USE ONLY

January 11, 2019