



ProSci Incorporated
12170 Flint Place
Poway, CA 92064, USA
prosci-inc.com

P: +1 (888) 513-9525
Local: +1 (858) 513-2638
Fax: +1 (858) 513-2692

customercare@prosci-inc.com
techsupport@prosci-inc.com

SARS-CoV-2 WT Pseudovirus


CATALOG NUMBER: 95-200

Specifications

Species	SARS-CoV-2 (COVID-19) Spike pseudovirus
Source Species	HEK293 cells
Tested Applications	Neut, Viral Infection

Properties

Physical State	Liquid
Buffer	DMEM/10% FBS
Storage Conditions	The pseudovirus is stable in a liquid state at -70° for 12 months. Avoid repeated freeze-thaw cycles.



SARS-CoV-2 WT Pseudovirus 1

Figure 1 Neutralization Assay of SARS-CoV-2 Pseudovirus by SARS-CoV-2 Spike Antibody
 Neutralization antibody: Anti-SARS-CoV-2 Spike S1 antibody, SD9853, diluti...

Disclaimer

Disclaimer	Products are intended for laboratory research purposes only and should be used by qualified personnel only. They are not intended for use in humans. ProSci is not liable for damages or injuries resulting from receipt and/or use of ProSci materials. Please refer to the Material Safety Data Sheet (MSDS) for safe storage, handling, and use procedures.
-------------------	--

For research use only. For additional information, visit ProSci's [Terms and Conditions Page](#).